



PurePedi

## USER MANUAL

Model: PECALREM



Please read this entire manual  
before using this product.


### INTRODUCTION

Thank you for purchasing PurePedi! This sleek, cordless callus solution is the perfect grooming tool to keep your feet soft and smooth to the touch.

### PRODUCT FEATURES

We've designed PurePedi to give you a safe yet effective alternative to harsh, metal scrapers. Here is an overview of the various features you will find on your PurePedi:

1. Roller with fine grain (Purple)
2. Roller holder
3. Slide switch (on/off switch)
4. Safety button
5. Cordless Callus Remover
6. Release button for the roller holder
7. Battery compartment lid
8. Roller with coarse grain (turquoise)
9. Protective lid
10. Cleaning Brush



### OPERATING YOUR PurePedi

#### Power Supply

PurePedi is powered by 2 AA batteries. To ensure the device is powered properly:

1. Slide the battery compartment lid to open.
2. Insert 2 AA batteries using polarity guides.
3. Slide the battery compartment lid to close.

#### Attaching and Removing Rollers

PurePedi comes with two different rollers. The turquoise roller has a coarse grain that removes thick calluses. The purple roller has a finer grain that removes calluses and smoothes skin.

Place the roller holder (labeled #2 in the diagram) on top of PurePedi (#5) so that the wheels of the roller and the PurePedi are aligned.

1. Press the roller holder on to PurePedi until you hear it connect.

2. Make sure it is secure.

3. To remove, press the release button (#6) and remove the roller holder using an upward motion.

#### Replacing Rollers

Your rollers should last several grooming sessions. However, you should immediately replace your rollers if you see signs of damage or it takes more time to achieve your normal results.

### REMOVING CALLUSES WITH THE PurePedi

1. Clean and dry the area you plan to treat.
2. Lay a blanket or towel down to catch removed calloused skin.
3. Press the safety button (#4 on the diagram) down and then press the slide switch (#3) up. The roller will begin to rotate.
4. Apply the roller slowly and carefully to the callus you wish to remove for 2-3 seconds.
5. Stop and check to see if you have achieved the desired effect.
6. If not, repeat step 4 until smooth.
7. Clean skin with warm water or a dry towel to remove excess skin.
8. Apply body lotion (optional).

**Note:** To avoid inflammation, never apply the PurePedi for more than 2-3 seconds at a time before stopping to check the results.

### CLEANING AND MAINTENANCE

#### After Each Use

1. Switch off the device.
2. Use the cleaning brush to clean the roller holder and the rollers.
3. Place protective lid on the roller to store.

#### From time to time

1. Wipe down the housing of the device with a damp cloth when necessary.
2. Check the surface of the roller for damage and replace if necessary.

### CLEANING AND MAINTENANCE (cont'd.)

DO NOT immerse any parts of PurePedi under water to clean.

DO NOT use astringent or abrasive cleaning agents.

### IMPORTANT HEALTH & SAFETY CONSIDERATIONS

This device should only be used for personal, home-based callus removal. It is not intended for commercial use.

DO NOT use on wet, reddened, or injured skin.

DO NOT use on healthy skin that is free from calluses.

Birth marks, warts and other protrusions should not be treated with PurePedi.

Consult with your physician prior to using the PurePedi if you have diabetes or high blood pressure.

DO NOT switch PurePedi on if the roller holder has been removed.

DO NOT allow children to play with PurePedi.

The cleaning brush is a small part that could be swallowed. Keep it away from children under the age of 3.

### WARRANTY INFORMATION

Pure Enrichment offers a 2-Year Warranty on your PurePedi. Visit [pureenrichment.com](http://pureenrichment.com) to contact us regarding warranty requests or questions.

### PRODUCT REGISTRATION

Registering your PurePedi ensures your 2-Year Warranty and allows you to receive product updates and streamlined customer support.

Register today at:  
[pureenrichment.com/productregistration](http://pureenrichment.com/productregistration)



Pure Enrichment™  
Santa Ana, CA 92705  
[pureenrichment.com](http://pureenrichment.com)