



**Osaki**



**Black**



**Tan**



**Brown**



**Charcoal**

The **All New OS-7200H** is the latest and most advanced massage chair in the Osaki series. Equipped with a **micro-computer**, it is engineered to perform various intelligent massage styles and combinations with precision. Accompanied with the new technology is the added pleasures of **heat therapy**.

## Computer Body Scan

Allows the rollers to scan the curvature of the user's body and self-adjust making each massage experience unique to your body type.



## Zero Gravity

Reclined position designed to relieve back pressure and enhance the massage experience.



## Quad Roller Head Massage System

Twin rollers to mimic the feel of human thumbs & fingers, enhancing the massage experience.



## 51 Air Bag Massage

Strategically placed bags inflate and massage vital areas such as: neck, shoulder, back, arms, hands, calves, feet, and the lumbar area.



## Auto Leg Scan

Auto adjusts to individuals' leg length to ensure the airbags are concentrated on the correct areas.



## Remote Control

Includes 6 unique preset programs and controls for manual control of methods, speed, intensity, and roller head width.



## Chromo-Therapy

Blue LED ambient lights create a chromo-therapeutic environment when in a dark room



## Multi-Layer Pillow & Padding / Removable Head Cover

Dual pad system allows the user to control the intensity of the massage: removing the pad stops the massage. The removable head cover can be hand washed.



## Hip & Pelvid Massage

Two sets of air bags placed on the sides of the seat compress the thighs and hips for an inwards massage.



## Outer Shoulder Massage

Located along the outer shoulders, air bags compress and apply inward pressure along the back of the neck and shoulders to sooth and loosen neck tension.



## Arm & Hand Massage

Multiple air bags are located on the top and bottom of the armrest that apply a compression massage along the forearm and hand.



## Calf & Foot Massage

Multiple air bags along the legs and feet inflate to apply a compression massage to the calves and feet. Pressure point balls along the bottom of the foot massager administer a concentrated pin-point massage.



## Stretch Mode

The air bags inflate to firmly grip the calves and slowly extend forward to stretch and pull down on the legs. This stretches the muscles and joints from the ankles to the lower back.



**Forward Stretch**  
The calves are held firmly and then extended out.



**Downward Stretch**  
The calves are held firmly and then pulled down.

## Auto Timer

Auto-Off timer can be set at 5-30 minutes and ensures a user that falls asleep during a session isn't affected by over use.

## iPod & MP3 Attachment

Play your favorite music from your device through the quality compact speakers mounted in the shoulder massagers.



## Heat Therapy

Located along the shoulder are heating pads that loosen muscle tissue: this relieves mild stress and tension and provides a deeper and more effective massage session.

