




## Grabcessories 2-in-1 Grab Bar/Toilet Paper Holder with Multiple EZ Clean Rubber Nubby Grips & FREE MOUNT DADDY Fasteners (2)


### Model #'s

- 61022 - Chrome
- 61028 - Brushed Nickel
- 61040 - Matte Black
- 2-in-1 design combines safety and functionality
- 1.25" diameter bar supports up to 500 lbs. in wall substrates that support the load
- Includes (2) FREE MOUNT DADDY Hollow Wall Fasteners
- Includes (2) LiveSafe Stud Mount Flanges designed to securely mount 3 screws per stud
- Stainless steel - no rust
- Multiple EZ Clean Rubber Nubby Grips on back of bar



### Installation Location Instructions

The 1st mounting post closest to the user when sitting on the toilet should be mounted at 41" from the back wall behind the toilet tank and 28.5" off the floor.



**NOTE:** This may look odd at first vs. normal paper holder placement but no worries as it provides residents optimal ability to pull forward elevating from seated position. If the bar is next to their side like a regular paper holder benefit is minimized